






Rehabilitation Aids & Orthopaedic Soft Goods



NOTES

|



CERVICAL SPLINTS

1-6



SHOULDER, RIB & CLAVICLE SPLINTS

7-12



BACK & POSTURE SPLINTS

13-20



ABDOMINAL SPLINTS

21-26



FINGER, WRIST, ELBOW & ARM SPLINTS

27-44



KNEE AND ANKLE SPLINTS

45-59



TRACTION KITS

60-64





ORTHOPAEDIC GENERAL AIDS & GEL PADS

65-77



CERVICAL SPLINTS






SIZE CHART

Getting the right splint starts with the right measurement. Our splint are custom made to fit a variety of clients.

- Using a tape measure where indicated arrow.
- Double check the accuracy of the measurement.

* Be sure limb is relaxed during the measurement.



HOW TO MEASURE:

-Measure the circumference around your neck & choose your desired neck collar height.

	XS	S	M	L	XL	XXL	U
Neck Circum. (Inch)	—	13-15	16-18	18-20	20-22	—	Universal
Neck Circum. (cm)	—	33-38	41-46	46-51	51-56	—	Universal
Height (Inch)	—	3	3.5	4	4.5	—	Universal
Height (cm)	—	6.5	7.6	10.2	11.4	—	Universal

* Customized sizes can be made as per requirement.

* This size chart is designed to be used as a general guide reference.

Cervical Collar (Soft) with Eyelet-Beige



FEATURES

- Unique contour of soft foam accommodates the chin comfortably.
- Double layered foam for warm comfort and snug fit around neck.
- Splint inside the collar provides stiffness for better neck support.
- Fine, knitted fabric cover absorbs sweat and is skin-friendly.
- Eyelets provide ventilation.

INDICATED FOR / APPLICATION

- Preventive care of neck during work, travel, watching TV or reading
- Minor neck injuries like strain, sprain and automobile whiplash
- Night-time use in combination with day-time use of semi-rigid or hard collar

Size	Item Code
S	RA5100S
M	RA5100M
L	RA5100L
XL	RA5100XL

Cervical Collar (Soft) with Eyelet-Grey



FEATURES

- Unique contour of soft foam accommodates the chin comfortably.
- Double layered foam for warm comfort and snug fit around neck.
- Splint inside the collar provides stiffness for better neck support.
- Fine, knitted fabric cover absorbs sweat and is skin-friendly.
- Eyelets provide ventilation.

INDICATED FOR / APPLICATION

- Preventive care of neck during work, travel, watching TV or reading
- Minor neck injuries like strain, sprain and automobile whiplash
- Night-time use in combination with day-time use of semi-rigid or hard collar

Size	Item Code
S	RA5101S
M	RA5101M
L	RA5101L
XL	RA5101XL

Cervical Collar Full Foam



FEATURES

- Soft foam collar with unique contour and fine resiliency offers warm comfort and snug fit around the neck.
- Fine, natural cotton stockinet cover absorbs sweat and remains skin friendly.

INDICATED FOR / APPLICATION

- Preventive care, to relieve muscle tension during work, travel, watching TV and reading
- Minor neck injuries like strain, sprain and automobile whiplash
- Night-time use along with daytime use of semi-rigid or hard collar

Size	Item Code
S	RA5105S
M	RA5105M
L	RA5105L
XL	RA5105XL

Cervical Collar (Hard)

FEATURES

- Made of high-density polythene sheet.
- Banian foam padding for comfort.
- Velcro closures to easy application removal and adjustment.
- Holes provided for air circulation.

INDICATED FOR / APPLICATION

- Cervical disc disease
- Cervical spondylosis
- Stabilization or hyperextension of the cervical region
- Torticollis

Size	Item Code
S	RA5110S
M	RA5110M
L	RA5110L
XL	RA5110XL



Cervical Collar Hard Adjustable

FEATURES

- Made of high-density polythene sheet.
- Two piece height adjustments to maintain desired flexion or extension.
- Banian foam padding for comfort.
- Velcro closures to easy application removal and adjustment.
- Holes provided for air circulation.

INDICATED FOR / APPLICATION

- Cervical disc disease
- Cervical spondylosis
- Stabilization or hyperextension of the cervical region
- Torticollis

Size	Item Code
S	RA5112S
M	RA5112M
L	RA5112L
XL	RA5112XL



Cervical Immobiliser (Philadelphia Collar)

FEATURES

- Two piece, pre-moulded design for perfect and comfortable fit.
- Mandibular (front) and occipital (back) splints limits flexion and extension movements of cervical spine.
- Made of unique non-absorbent sturdy foam that is easy to clean and maintain.
- Frontal opening for tracheotomy.
- Ventilation holes add comfort for extended use.

INDICATED FOR / APPLICATION

- Spondylosis
- Hyper-extension injury
- Post-surgical care
- Post-halo rehabilitative support

Size	Item Code
S	RA5120S
M	RA5120M
L	RA5120L
XL	RA5120XL



Cervical Splints

Cervical Immobiliser (Sprain Collar)



FEATURES

- Single piece collar with pre-formed chin-up to control flexion.
- Sturdy non-absorbent foam construction.
- Large frontal opening allows tracheotomy.
- Ventilation holes keep patient neck cool and dry.

INDICATED FOR / APPLICATION

- Posture maintenance during activities of daily living
- Whiplash injuries

Size	Item Code
S	RA5122S
M	RA5122M
L	RA5122L

Ambulance Collar



FEATURES

- Single piece flat design.
- Semi_rigid plastic, padded with sturdy foam provides comfortable immobilization.

INDICATED FOR / APPLICATION

- It is used for temporary cervical immobilization for accidental victims during the transportation & diagnostic studies.
- It is used as temporary support to neck in accidental conditions. it provides rigid support to limit the neck movement or to prevent harmful neck movements.

Size	Item Code
S	RA5124S
M	RA5124M
L	RA5124L
XL	RA5124XL

Cervical Pillow



FEATURES

- Soft foam collar with unique contour and fine resiliency offers warm comfort and snug fit around the neck.
- Fine, natural cotton stockinet cover absorbs sweat and remains skin friendly.

INDICATED FOR / APPLICATION

- Preventive care, to relieve muscle tension during work, travel, watching TV and reading
- Minor neck injuries like strain, sprain and automobile whiplash
- Night-time use along with daytime use of semi-rigid or hard collar



Size	Item Code
S	RA5130S
M	RA5130M
L	RA5130L
XL	RA5130XL



FEATURES

- Soft foam collar with unique contour and fine resiliency offers warm comfort and snug fit around the neck.
- Fine, natural cotton stockinet cover absorbs sweat and remains skin friendly.

INDICATED FOR / APPLICATION



- Preventive care, to relieve muscle tension during work, travel, watching TV and reading
- Minor neck injuries like strain, sprain and automobile whiplash
- Night-time use along with daytime use of semi-rigid or hard collar

Size	Item Code
S	RA5131S
M	RA5131M
L	RA5131L
XL	RA5131XL





SHOULDER, RIB & CLAVICLE SPLINTS






Size Chart

Getting the right splint starts with the right measurement. Our splint are custom made to fit a variety of clients.

- Using a tape measure where indicated arrow.
- Double check the accuracy of the measurement.

* Be sure limb is relaxed during the measurement.



HOW TO MEASURE:

Find the circumference around your chest below underarms.

	XS	S	M	L	XL	XXL	XXXL
For Chest length (inch)	24-27	28-31	32-35	36-39	40-43	44-47	48-51
For Chest length (cm)	61-69	71-79	81-89	91-99	101-109	111-119	121-129

- Tolerance ± 1 Inch
- Circumference measurement taken around Chest below underarms.

* Customized sizes can be made as per requirement.

* This size chart is designed to be used as a general guide reference.

SHOULDER, RIB & CLAVICLE SPLINTS

Universal Shoulder Immobilizer (Deluxe)-Beige



FEATURES

- Made from durable material.
- Straps fit over shoulders, cross at the back and knot around waist for even weight distribution and maximum patient comfort.

INDICATED FOR / APPLICATION

- Recuperation from broken bone
- Sprains
- Surgery

Size	Item Code
S	RA5200S
M	RA5200M
L	RA5200L
XL	RA5200XL

Universal Shoulder Immobilizer (Deluxe)-Grey



FEATURES

- Made from durable material.
- Straps fit over shoulders, cross at the back and knot around waist for even weight distribution and maximum patient comfort.

INDICATED FOR / APPLICATION

- Recuperation from broken bone
- Sprains
- Surgery

Size	Item Code
S	RA5201S
M	RA5201M
L	RA5201L
XL	RA5201XL

Meek Clavicle Brace-Deluxe (Velcro Foam Padding)



FEATURES

- Heavy duty straps with fabric-covered foam pads form "figure of eight" compression over the clavicular region.
- Inter-scapulae pad made of fabric-covered foam acts as base for straps & plastic buckles.
- Simple double-buckle system allows easy strap length adjustment and firm locking.

INDICATED FOR / APPLICATION

- Acromio-clavicular and scapular separations due to sprain
- Management of clavicle injuries
- Posture correction in the case of stooped shoulders.

Size	Item Code
S	RA5210S
M	RA5210M
L	RA5210L
XL	RA5210XL

Meek Clavicle BraceDeluxe-Beige

FEATURES

- Heavy duty straps with fabric-covered foam pads form "figure of eight" compression over the clavicular region.
- Inter-scapulae pad made of fabric-covered foam acts as base for straps & plastic buckles.
- Simple double-buckle system allows easy strap length adjustment and firm locking.

INDICATED FOR / APPLICATION

- Acromio-clavicular and scapular separations due to sprain
- Management of clavicle injuries
- Posture correction in the case of stooped shoulders.

Size	Item Code
S	RA5212S
M	RA5212M
L	RA5212L
XL	RA5212XL



Meek Clavicle BraceDeluxe-Grey

FEATURES

- Heavy duty straps with fabric-covered foam pads form "figure of eight" compression over the clavicular region.
- Inter-scapulae pad made of fabric-covered foam acts as base for straps & plastic buckles.
- Simple double-buckle system allows easy strap length adjustment and firm locking.

INDICATED FOR / APPLICATION

- Acromio-clavicular and scapular separations due to sprain
- Management of clavicle injuries
- Posture correction in the case of stooped shoulders.

Size	Item Code
S	RA5213S
M	RA5213M
L	RA5213L
XL	RA5213XL



Universal Shoulder Immobilizer (Ring Type)-Beige

FEATURES

- High quality elastic for biceps & forearm cuffs provides support to shoulder without restricting respiratory chest movement.
- The biceps and forearm cuffs keep forearm across chest, limiting shoulder rotation & abduction in a comfortable way.
- Designed to fit either left or right arm.

INDICATED FOR / APPLICATION

- Soft tissue injuries like strain and sprain in the shoulder joint
- Shoulder dislocation
- Bursitis
- Emergency splinting to secure the arm & shoulder

Size	Item Code
S	RA5220S
M	RA5220M
L	RA5220L
XL	RA5220XL



SHOULDER, RIB & CLAVICLE SPLINTS

Universal Shoulder Immobilizer (Ring Type)-Grey



FEATURES

- High quality elastic for biceps & forearm cuffs provides support to shoulder without restricting respiratory chest movement.
- The biceps and forearm cuffs keep forearm across chest, limiting shoulder rotation & abduction in a comfortable way.
- Designed to fit either left or right arm.

INDICATED FOR / APPLICATION

- Soft tissue injuries like strain and sprain in the shoulder joint
- Shoulder dislocation
- Bursitis
- Emergency splinting to secure the arm & shoulder

Size	Item Code
S	RA5221S
M	RA5221M
L	RA5221L
XL	RA5221XL

Arm Sling Strap-Beige



FEATURES

- For holding injured arm, shoulder and wrist Heavy duty straps assures sturdy design Tubular neck pad increase comfort level Slide buckle provides maximum adjustability
- It is used in immobilization of hand, wrist & shoulder

Size	Item Code
S	RA5230S
M	RA5230M
L	RA5230L
XL	RA5230XL

Arm Sling Strap-Grey



FEATURES

- For holding injured arm, shoulder and wrist Heavy duty straps assures sturdy design Tubular neck pad increase comfort level Slide buckle provides maximum adjustability
- It is used in immobilization of hand, wrist & shoulder

Size	Item Code
S	RA5231S
M	RA5231M
L	RA5231L
XL	RA5231XL

Rib Belt (Female)-Beige

FEATURES

- Superior quality, breathable elastic band offers compression, support while allowing normal respiration.
- The expansion of the elastic is proportional to the expansion of the rib cage.
- Velcro closure enables easy application in proper alignment.

INDICATED FOR / APPLICATION

- Compression support
- To protect surgical incisions on chest
- To stabilize sternum fracture
- To limit painful chest expansion in the case of rib fracture and/or dislocation

Size	Item Code
S	RA5240S
M	RA5240M
L	RA5240L
XL	RA5240XL

RA 5242 Rib Belt (Male) Also available - Beige



Rib Belt (Female)-Grey

FEATURES

- Superior quality, breathable elastic band offers compression, support while allowing normal respiration.
- The expansion of the elastic is proportional to the expansion of the rib cage.
- Velcro closure enables easy application in proper alignment.

INDICATED FOR / APPLICATION

- Compression support
- To protect surgical incisions on chest
- To stabilize sternum fracture
- To limit painful chest expansion in the case of rib fracture and/or dislocation

Size	Item Code
S	RA5241S
M	RA5241M
L	RA5241L
XL	RA5241XL

RA 5243 Rib Belt (Male) Also available - Grey



Hyper Extension Brace (Ash Brace)

FEATURES

- Designed to act as three-point principle by exerting pressure against sternum and pubis on anterior aspects and the Lumbar region on the Posterior Aspect.
- This fixes the Lumbar vertebral column in a hyperlordotic position.
- Can be fitted to any patient due to the sliding arrangements.



INDICATED FOR / APPLICATION

- Diskettes
- Tuberculosis
- Undisplaced stable fractures
- Adolescent kyphosis
- Any problem of the thoracic Lumbar junction with needs moderate immobilization

Size	Item Code
Universal	RA 5250



BACK & POSTURE SPLINTS



SIZE CHART


Comfo Smart



Getting the right splint starts with the right measurement. Our splint are custom made to fit a variety of clients.

- Using a tape measure where indicated arrow.
- Double check the accuracy of the measurement.

* Be sure limb is relaxed during the measurement.



HOW TO MEASURE:

To take accurate measurement.

- Measure around your natural waist line, at the height of your navel
- Measure around your fullest part of your chest, keeping the tape parallel to the floor.
- Measure torso length starting from the back of neck(base) measure down to your natural waist.

Make sure you don't draw the tape too tightly.

Waist

	S	M	L	XL	XXL
For waist length (inch)	28-32	32-36	36-40	40-44	44-48
For waist length (cm)	71-81	81-91	91-101	101-111	111-121

NOTE: Sizes above are not jean(waist sizes)

Chest

For Chest length (inch)	28-31	32-35	36-39	40-43	44-47
For Chest length (cm)	71-79	81-89	91-99	101-109	111-119

Torso Height

For Torso length (inch)	21.5-22.5	22.5-23.5	23.5-24.5	24.5-25.5	25.5-26.5
For Torso length (cm)	55-57	57-60	60-63	63-66	66-69

* Customized sizes can be made as per requirement.

* This size chart is designed to be used as a general guide reference.

Cervical Splints

Shoulder, Rib & Clavicle Splints

Back & Posture Splints

Abdominal Splints

Finger, Wrist, Elbow & Arm Splints

Knee & Ankle Splints

Traction Kits

Orthopaedic General Aids & Gel Pads

Sacro Lumber Brace (Belt) Contoured-Beige



FEATURES

- Double elastic straps and thin material provides uniform compression and comfortable support to back section
- Four heavy aluminium malleable metal stays are removable & can be given desired contoured to support the specific vertebrae.
- Easy to apply, and easy to adjust different sizes as per user's requirements.
- It is beneficial for patients having pendulous abdomen with lower back problems, degenerative disc problems, lumbosacral spondylolysis, osteoporotic pains of lumbar spines.
- To provide uniform support to lower back section with comfortable immobilization.

Size	Item Code
S	RA5300S
M	RA5300M
L	RA5300L
XL	RA5300XL

Sacro Lumber Brace (Belt) Contoured-Grey



FEATURES

- Double elastic straps and thin material provides uniform compression and comfortable support to back section
- Four heavy aluminium malleable metal stays are removable & can be given desired contoured to support the specific vertebrae.
- Easy to apply, and easy to adjust different sizes as per user's requirements.
- It is beneficial for patients having pendulous abdomen with lower back problems, degenerative disc problems, lumbosacral spondylolysis, osteoporotic pains of lumbar spines.
- To provide uniform support to lower back section with comfortable immobilization.

Size	Item Code
S	RA5301S
M	RA5301M
L	RA5301L
XL	RA5301XL

Sacro Lumber Brace (Belt) Contoured Economy-Beige



FEATURES

- Double elastic straps and thin material provides uniform compression and comfortable support to back section
- Four thin aluminium malleable metal stays are removable & gives effective support to the back.
- Easy to apply, and easy to adjust different sizes as per user's requirements.
- It is beneficial for patients having pendulous abdomen with lower back problems, degenerative disc problems, lumbosacral spondylolysis, osteoporotic pains of lumbar spines.
- To provide uniform support to lower back section with comfortable immobilization

Size	Item Code
S	RA5302S
M	RA5302M
L	RA5302L
XL	RA5302XL

Sacro Lumber Brace (Belt) Contoured Economy-Grey

FEATURES

- Double elastic straps and thin material provides uniform compression and comfortable support to back section
- Four thin aluminium malleable metal stays are removable & gives effective support to the back.
- Easy to apply, and easy to adjust different sizes as per user's requirements.
- It is beneficial for patients having pendulous abdomen with lower back problems, degenerative disc problems, lumbosacral spondylolysis, osteoporotic pains of lumbar spines.
- To provide uniform support to lower back section with comfortable immobilization

Size	Item Code
S	RA5303S
M	RA5303M
L	RA5303L
XL	RA5303XL



Sacro Lumber Brace (Belt) Super Fine

FEATURES

- Fully elasticated to add extra pressure and support to the sacro & lumbar regions.
- Splints made of tampered steel strips / plastic strips.
- Easy to apply, and easy to adjust different sizes as per user's requirements.
- It is beneficial for patients with lower back problems, degenerative disc problems, lumbosacral spondylolysis, osteoporotic pains of lumbar spines.
- To provide uniform support to lower back section with rigid immobilization.

Size	Item Code
S	RA5310S
M	RA5310M
L	RA5310L
XL	RA5310XL



Sacro Lumber Brace (Belt) Super Fine Double lock

FEATURES

- Fully elasticated to add extra pressure and support to the sacro & lumbar regions.
- Double elastic support for extra support.
- Splints made of tampered steel strips / plastic strips.
- Easy to apply, and easy to adjust different sizes as per user's requirements.
- It is beneficial for patients with lower back problems, degenerative disc problems, lumbosacral spondylolysis, osteoporotic pains of lumbar spines.
- To provide uniform support to lower back section with rigid immobilization.

Size	Item Code
S	RA5312S
M	RA5312M
L	RA5312L
XL	RA5312XL



Sacro Lumber Brace (Belt) Fine



FEATURES

- Fine quality elastic to add extra pressure and support to the sacro & lumbar regions.
- Splints made of tampered steel strips / plastic strips.
- Easy to apply, and easy to adjust different sizes as per user's requirements.
- It is beneficial for patients with lower back problems, degenerative disc problems, lumbosacral spondylolysis, osteoporotic pains of lumbar spines
- To provide uniform support to lower back section with rigid immobilization.

Size	Item Code
S	RA5320S
M	RA5320M
L	RA5320L
XL	RA5320XL

Sacro Lumber Brace (Belt) Fine Double lock



FEATURES

- Fine quality elastic to add extra pressure and support to the sacro & lumbar regions.
- Double elastic support for extra support.
- Splints made of tampered steel strips / plastic strips.
- Easy to apply, and easy to adjust different sizes as per user's requirements.
- It is beneficial for patients with lower back problems, degenerative disc problems, lumbosacral spondylolysis, osteoporotic pains of lumbar spines
- To provide uniform support to lower back section with rigid immobilization.

Size	Item Code
S	RA5322S
M	RA5322M
L	RA5322L
XL	RA5322XL

Eco Frame Back Support



FEATURES

- Light weight metal frame with flexibility to adjust the contoured of the back region
- Special elastic and thin material provides uniform compression and scientific design gives rigid support to back section
- Easy to apply, and easy to remove. Due to hook and loop arrangement
- Complete round back support of metal frame restricts the side movement of the back
- It is useful in treatment of acute lumbago, slip disc recurrent lower back
- To provide firm support to back and lumbar region

Size	Item Code
S	RA5330S
M	RA5330M
L	RA5330L
XL	RA5330XL

Sacro-Lumber Belt (Contoured Type)

FEATURES

- Constructed of superior quality perforated elastic for extra breathability.
- Back panel with heavy duty, durable fabric and soft pad for support to lumbo-sacral region.
- Four flexible stays provide rigid support and prevent edge roll over.
- Elastic front panel in tapered style ensures proper fit in hip/waist region.
- Designed to encourage good body mechanism.

INDICATED FOR / APPLICATION

- Ideal as a safety support against work-related back injuries especially for those handling heavy weights

Size	Item Code
S	RA5340S
M	RA5340M
L	RA5340L
XL	RA5340XL



Lumbo sacral support Premium-Beige

FEATURES

- Lumbo sacral support is designed to provide comfortable support to the lumbar and sacral region of back. It helps in mild back ache.
- Breathable superior quality porous elastic offers good compression support. Flexible splints on back for proper support to the back region. Wide hook and loop closure for easy application. Double pull mechanism for better fitting around the waist.
- It is used in support to back section for mild back ache and for post operative care. It is ideal for everyday use to support at home or while working in office.

Size	Item Code
S	RA5350S
M	RA5350M
L	RA5350L
XL	RA5350XL



Lumbo sacral support Premium-Grey

FEATURES

- Lumbo sacral support is designed to provide comfortable support to the lumbar and sacral region of back. It helps in mild back ache.
- Breathable superior quality porous elastic offers good compression support. Flexible splints on back for proper support to the back region. Wide hook and loop closure for easy application. Double pull mechanism for better fitting around the waist.
- It is used in support to back section for mild back ache and for post operative care. It is ideal for everyday use to support at home or while working in office.

Size	Item Code
S	RA5351S
M	RA5351M
L	RA5351L
XL	RA5351XL



BACK & POSTURE SPLINTS

Sternal Splint



FEATURES

- Special soft foam pad provides mild compression support to sternum
- Easy to apply and remove
- Allows normal respiration
- It is useful in protecting surgical incision on chest. for protection of sternum after heart surgery

Size	Item Code
S	RA5360S
M	RA5360M
L	RA5360L
XL	RA5360XL

Posture Aid



FEATURES

- A comfortable aid for patients with stooping shoulders.
- Holds shoulder upright and keeps back in correct posture with chest out position.

Size	Item Code
S	RA5370S
M	RA5370M
L	RA5370L
XL	RA5370XL

Posture aid magnetic support



Features

- Supplying gentle back support.
- Perfect posture correcting device.
- Available for both men and women.
- Pulling shoulder and back to correct posture.
- Fully adjustable, comfortable and easy to wear.
- High energetic magnets contact the spine and lumbar region.
- Designed to help correct wrong body posture of back and shoulder.
- Capable of preventing customary humpback adolescent, giving teenage a health and beautiful youth.

Size	Item Code
S	RA5372S
M	RA5372M
L	RA5372L
XL	RA5372XL

Taylor's Brace Dorso Lumbar Spinal Brace (Super Deluxe)-Beige

FEATURES

- Two rigid upright para-spinal and lateral bars of light-weight are pre-contoured to keep spine in it's natural position.
- Elastic abdominal panels provide the required compression.
- Thoracic straps of 4" elastic provide comfortable fixation of the brace to the body.
- Clavicle straps with cotton stockinet covered foam pads correct shoulder posture.

INDICATED FOR / APPLICATION

- Mild thoraco-lumbar injuries • Herniated disc problems • Post-operative rehabilitation

Size	Item Code
S	RA5380S
M	RA5380M
L	RA5380L
XL	RA5380XL



Taylor's Brace Dorso Lumbar Spinal Brace (Super Deluxe)-Grey

FEATURES

- Two rigid upright para-spinal and lateral bars of light-weight are pre-contoured to keep spine in it's natural position.
- Elastic abdominal panels provide the required compression.
- Thoracic straps of 4" elastic provide comfortable fixation of the brace to the body.
- Clavicle straps with cotton stockinet covered foam pads correct shoulder posture.





INDICATED FOR / APPLICATION

- Mild thoraco-lumbar injuries • Herniated disc problems • Post-operative rehabilitation

Size	Item Code
S	RA5381S
M	RA5381M
L	RA5381L
XL	RA5381XL



ABDOMINAL SPLINTS




SIZE CHART

Getting the right splint starts with the right measurement. Our splint are custom made to fit a variety of clients.

- Using a tape measure where indicated arrow.
- Double check the accuracy of the measurement.

* Be sure limb is relaxed during the measurement.



HOW TO MEASURE:

Measure around your natural waist line, at the height of your navel

S M L XL XXL

For waist length
(inch)

28-32 32-36 36-40 40-44 44-48

For waist length
(cm)

71-81 81-91 91-101 101-111 111-121

* Customized sizes can be made as per requirement.

* This size chart is designed to be used as a general guide reference.

Cervical Splints

Shoulder, Rib &
Clavicle Splints

Back & Posture
Splints

Abdominal
Splints

Finger, Wrist, Elbow
& Arm Splints

Knee & Ankle
Splints

Traction Kits

Orthopaedic General
Aids & Gel Pads

ABDOMINAL SPLINTS

6" Adbominal Support



FEATURES

- Superior quality, breathable elastic band of 10" width offers good compression support.
- Flexible stays prevent edge roll-over on usage.
- Velcro closure for easy application.

INDICATED FOR / APPLICATION

- Post-surgical compression support in abdominal region

Size	Item Code
S	RA5400S
M	RA5400M
L	RA5400L
XL	RA5400XL

8" Abdominal Support - Premium-Beige



FEATURES

- Abdominal Support / Tummy trimmer Premium is designed for Abdomen support and compresses the muscle for waist line reduction.
- Three panel superior porous elastic improves ventilation and patient comfort.
- Broad hook loop panel allows controlled compression.
- It is better for sizing flexibility.
- Easy to wear
- It is beneficial for patients with back problems, post operative supports to abdominal region, for mothers after delivery.

Size	Item Code
S	RA5410S
M	RA5410M
L	RA5410L
XL	RA5410XL

8" Abdominal Support - Premium-Grey



FEATURES

- Abdominal Support / Tummy trimmer Premium is designed for Abdomen support and compresses the muscle for waist line reduction.
- Three panel superior porous elastic improves ventilation and patient comfort.
- Broad hook loop panel allows controlled compression.
- It is better for sizing flexibility.
- Easy to wear
- It is beneficial for patients with back problems, post operative supports to abdominal region, for mothers after delivery.

Size	Item Code
S	RA5411S
M	RA5411M
L	RA5411L
XL	RA5411XL

9" Abdominal Support

FEATURES

- Abdominal Support 9 inches single elastic is designed for abdomen support and compress the muscles for waist line reduction and comfortable support to back section.
- 9 inches single elastic improves ventilation
- Strong and porous elastic webbing ensures good compression
- No roll over, easy to apply
- Broad hook and loop panel allowed controlled compression
- Nine inches single elastic offers good compression support, improves ventilation and comfort. Flexible stays prevent roll over on usage. Wide hook and loop closure for easy application.
- It is used in compression to abdomen muscles and support to back section. It is used in post operative and post delivery aids

Size	Item Code
S	RA5412S
M	RA5412M
L	RA5412L
XL	RA5412XL



9" Abdominal Support

FEATURES

- Abdominal Support 9 inches single elastic is designed for abdomen support and compress the muscles for waist line reduction and comfortable support to back section.
- 9 inches single elastic improves ventilation
- Strong and porous elastic webbing ensures good compression
- No roll over, easy to apply
- Broad hook and loop panel allowed controlled compression
- Nine inches single elastic offers good compression support, improves ventilation and comfort. Flexible stays prevent roll over on usage. Wide hook and loop closure for easy application.
- It is used in compression to abdomen muscles and support to back section. It is used in post operative and post delivery aids

Size	Item Code
S	RA5413S
M	RA5413M
L	RA5413L
XL	RA5413XL



8" Abdominal Support Fine (Full Elastic) Tummy Trim Belt

FEATURES

- Fine quality, breathable elastic band offers good compression support.
- Flexible stays prevent edge roll-over on usage.
- Velcro closure for easy application.

INDICATED FOR / APPLICATION

- Post-surgical compression support in abdominal region

Size	Item Code
S	RA5420S
M	RA5420M
L	RA5420L
XL	RA5420XL



8" Abdominal Support Super Fine (Full Elastic) Tummy Trim Belt



FEATURES

- Superior quality, breathable elastic band offers good compression support.
- Flexible stays prevent edge roll-over on usage.
- Velcro closure for easy application.

INDICATED FOR / APPLICATION

- Post-surgical compression support in abdominal region

Size	Item Code
S	RA5422S
M	RA5422M
L	RA5422L
XL	RA5422XL

8" Abdominal Support Neo Type Black



FEATURES

- The Neo type- drytex abdominal support is constructed using a unique polyester nylon cotton weave three layered fabric that provides warmth compression and support to the abdomen. It is used in compression to abdomen muscles and support to back section
- Allows for a wider range of adjustability and easy application.
- Sizes : UNIVERSAL (28" - 44" length) / UNIVERSAL (28" – 54" length)

Size	Item Code
Universal	RA 5430

9" Abdominal Support (Double Lock)



FEATURES

- Superior quality, breathable elastic band offers good compression support.
- Flexible stays prevent edge roll-over on usage.
- Velcro closure for easy application.

INDICATED FOR / APPLICATION

- Post-surgical compression support in abdominal region

Size	Item Code
S	RA5440S
M	RA5440M
L	RA5440L
XL	RA5440XL

10" Abdominal Support

FEATURES

- Superior quality, breathable elastic band of 10" width offers good compression support.
- Flexible stays prevent edge roll-over on usage.
- Velcro closure for easy application.

INDICATED FOR / APPLICATION

- Post-surgical compression support in abdominal region

Size	Item Code
S	RA5451S
M	RA5451M
L	RA5451L
XL	RA5451XL



Maternity support

Applications

- Provides support to the lumbar and pelvic region.
- Relieves the strain on muscles, ligaments caused by the pendulous abdomen.
- Helps with proper posture and relieves stress to the baby by encouraging proper body position and movement.
- Soft pile bonded to foam
- Helps relieve pain by gently lifting the abdomen.

Size	Item Code
Universal	RA 5460



Maternity Belt with Back Support

Features:

- It is an effective support system that effectively helps to reduce lower back pain and fatigue during pregnancy.
- The Brace provides effective support to lower lumbar region, abdomen, and hip area to provide relief from pain and discomfort during pregnancy.

- Fully adjustable to accommodate abdominal growth throughout pregnancy.
- Constructed of a multi-layered laminate with a cotton/nylon lining.
- Four elastic side panels add compression and ventilation.
- Included abdominal lift attachment gives lift without unnecessary pressure.
- Built-in support holds the belt in place and prevents the belt from rolling.

Size	Item Code
Universal	RA 5462





FINGER, WRIST, ELBOW & ARM SPLINTS




SIZE CHART

Getting the right splint starts with the right measurement. Our splint are custom made to fit a variety of clients.

- Using a tape measure where indicated arrow.
- Double check the accuracy of the measurement.

* Be sure limb is relaxed during the measurement.



ELBOW

HOW TO MEASURE:

- Measure the circumference at elbow joint.
- Measure the circumference at the wrist joint.



WRIST

	XS	S	M	L	XL	XXL	U
WRIST							
Wrist Circum. (Inch)	—	6"	6.5"	7"	7.5"	8"	Universal
Wrist Circum. (cm)	—	15	16.5	17.5	19	20	Universal
ELBOW							
Elbow Circum. (Inch)	—	8-10"	10.5-12"	12.5-14"	14.5-16"	16.5-17"	Universal
Elbow Circum. (cm)	—	20-25	26.5-30	31.5-35.5	36.5-40.5	41.5-43	Universal

* Customized sizes can be made as per requirement.

* This size chart is designed to be used as a general guide reference.

Cervical Splints

Shoulder, Rib & Clavicle Splints

Back & Posture Splints

Abdominal Splints

Finger, Wrist, Elbow & Arm Splints

Knee & Ankle Splints

Traction Kits

Orthopaedic General Aids & Gel Pads

FINGER, WRIST, ELBOW & ARM SPLINTS

Dynamic Finger Flexion Splint

FEATURES

Springs attached to give extension.

For extra comfort proximal and distal pads are covered with soft padding which ensure minimal finger slippage.

Flexes the pipjoint without interfering the mcp joint



Size	Item Code
S	RA5500S
M	RA5500M
L	RA5500L

Finger Cot

FEATURES

- Plastic coated malleable
- Aluminium can be moulded to keep the interphalangeal joints in required position.
- Poly foam padding for comfort.
- Fits any finger.
- Can be molded by bare hands.

INDICATED FOR / APPLICATION

- Immobilizes the fingers following reduction of fractures of proximal interphalanges • In extensor tendon injuries of the fingers • Collateral ligament injuries of the interphalangeal joint • Boutonniere's injury



Size	Item Code
S	RA5502S
M	RA5502M
L	RA5502L

Finger cot with velcro

FEATURES

Design allows application without taping

U shaped alluminium splint covers finger tip and assures rigid immobilization

Special foam padding provides comfortable fitting

used for fingure tip injuries and burn injuries

To protect finger in finger tip injuries,nail bed injuries



Size	Item Code
S	RA5504S
M	RA5504M
L	RA5504L

Finger Protector Splint with velcro

FEATURES

- Provide complete immobilization & protection for the finger.
- Ideal for treatment of finger fracture, sprains, dislocations or as protection following injury.
- Provide more comfort & accurate fitting.

Size	Item Code
Universal	RA 5506



Finger Protector Splint

FEATURES

To protect injured finger tip
Can be bent at any degree
Easy to apply and remove
It is useful in finger tip injuries

Size	Item Code
Universal	RA 5508



Straight Splint (Mallet Splint)

FEATURES

- Simple flat design to protect finger.
- Size allows use for full length or part of the finger.

INDICATED FOR / APPLICATION

• Used in the treatment of degenerative arthritis, pulp fractures and soft tissue injuries • Contractures due to burn injuries • Conditions requiring support and immobilization of interphalangeal joints • Dislocations or sprain during sports and crushed fractures due to accidents

Size	Item Code
S	RA5510S
M	RA5510M
L	RA5510L



FINGER, WRIST, ELBOW & ARM SPLINTS

Finger Mallet Splint



FEATURES

- For correction of mallet finger deformity
- Easily malleable to get desired shape
- Provides strong immobilization
- Easy to apply and remove
- It is useful in treatment of mallet finger deformity.

Size	Item Code
Universal	RA5512

Stax Mallet



FEATURES:

- Made from engineering plastic that is durable & contains ventilations holes.
- Splints keeps the finger in 10 degree of hyper extension
- It is used in hyper flexion injuries like mallet finger

Size	Item Code
Universal	RA5514

Frog Splint



FEATURES

- Plastic coated malleable.
- Aluminium of ease of forming and protective foam padding.
- Fits any finger.
- Can be molded by bare hands.

INDICATED FOR / APPLICATION

- Mallet finger • To hold interphalangeal joint in correct alignment following distal joint injuries • To keep distal interphalangeal joint in hyperextension

Size	Item Code
S	RA5516S
M	RA5516M
L	RA5516L

Finger Baseball Splint

FEATURES

- Innovative finger tip fold over design
- Alluminium splint can be reshaped as per required finger position
- Special foam padding provides comfortable fitting.
- Hyper flexion finger injuries due to the direct hit of the balls. in the game like cricket, baseball etc

Size	Item Code
S	RA5518S
M	RA5518M
L	RA5518L



Finger Extension Splint

FEATURES

- Plastic coated malleable
- Aluminium can be molded to keep the interphalangeal joints in required position.
- Poly foam padding for comfort.
- Fits any finger.
- Can be molded by bare hands.

INDICATED FOR / APPLICATION

- Immobilizes the fingers following reduction of fractures of proximal interphalanges • In extensor tendon injuries of the fingers • Collateral ligament injuries of the interphalangeal joint • Boutonniere's injury

Size	Item Code
S	RA5520S
M	RA5520M
L	RA5520L



Finger Extension Splint with Velcro

FEATURES

- Light weight splints provides effective immobilization of the fingers
- Can be moulded repeatedly.
- Foam padding provides comfort & ventilation to the finger
- easy to apply by just inserting the finger.
- Immobilization of the fingers following fracture. Direct hit of the balls. in the game like cricket, baseball etc

Size	Item Code
S	RA5522S
M	RA5522M
L	RA5522L



FINGER, WRIST, ELBOW & ARM SPLINTS

Thumb Abduction Splint



FEATURES

- For sprains & strains to the thumb's metacarpal joint.
- Special design minimize thumb movements
- Easy to apply and remove
- Useful for skier's thumb, bowler's thumb and general thumb injuries

Size	Item Code
S	RA5524S
M	RA5524M
L	RA5524L

Knuckle Binder



FEATURES

- Concepts of circular stainless steel spring are applied to maintain constant pressure unlike rubber band.
- It is well padded, unbreakable, washable, reusable and economic.

INDICATED FOR / APPLICATION

- To flex the metacarpophalangeal joints of the hand
- The stainless steel springs are detachable to enable the increase and decrease in tension of the apparatus as per the patient's requirement

Size	Item Code
S	RA5526S
M	RA5526M
L	RA5526L

Carpal tunnel splint



FEATURES

- To prevention of ct syndrome occurrence.
- Elastic band with velcro closure keeps palmar splint in position & gives required compression.
- Thumb loop allows self application of the product.
- High quality felt pad with rigid but light metal splint act as palmar support to keep the wrist in natural position.

Size	Item Code
Universal	RA 5528

Thumb Spica Splint

FEATURES

- Full hand and finger movement is possible
- Reduces pain and inflammation
- Customised abduction of thumb is possible
- For thumb injuries
- To support the cmc joint and immobilize the mcp joint of thumb without affecting hand movement

AVAILABLE SIZES

Universal

Size	Item Code
S	RA5530S
M	RA5530M
L	RA5530L



Elastic Wrist Brace-Beige

FEATURES

- Simple to apply, easy to adjust.
- Just touch to close.
- Peel to open.
- High durability elastic material.

INDICATED FOR / APPLICATION

- For strains and strains post fracture, bursitis microtrauma
- Tenosynovitis, arthritis and other conditions • Protects both wrist and upper bone

Size	Item Code
Universal	RA 5532



Elastic Wrist Brace-Grey

FEATURES

- Simple to apply, easy to adjust.
- Just touch to close.
- Peel to open.
- High durability elastic material.

INDICATED FOR / APPLICATION

- For strains and strains post fracture, bursitis microtrauma
- Tenosynovitis, arthritis and other conditions • Protects both wrist and upper bone

Size	Item Code
Universal	RA 5533



FINGER, WRIST, ELBOW & ARM SPLINTS

Wrist Brace (Neo type)



Wrist brace supports, protects and immobilizes the wrist while allowing free and natural movement of the thumb and fingers.

- Versatile Velcro closure.
- Controlled four way compressions.
- Free thumb and finger movement.

FEATURES

It is made of 3 layered fabric, which is four way stretchable. It provides more comfort, good grip and better durability. Wrap around design and fabric offers a versatile closing system, ensures easy wearing and removal. It allows customized compression, enhanced grip, better sizing flexibility, good fitting, better support and wrist stabilization.

Size	Item Code
Universal	RA 5534

Wrist Binder / Wrap-Beige



FEATURES

- Breathable, high quality elastic band offers mild compression over the wrist joint.
- Velcro closure enables easy application

INDICATED FOR / APPLICATION

- Non-surgical care of strained or sprained wrist joint • Preventive care for the wrist during activities like sports, driving, household work and assembling jobs

Size	Item Code
S	RA5536S
M	RA5536M
L	RA5536L

Wrist Binder / Wrap-Grey



FEATURES

- Breathable, high quality elastic band offers mild compression over the wrist joint.
- Velcro closure enables easy application

INDICATED FOR / APPLICATION

- Non-surgical care of strained or sprained wrist joint • Preventive care for the wrist during activities like sports, driving, household work and assembling jobs

Size	Item Code
S	RA5537S
M	RA5537M
L	RA5537L

Wrist Binder With Double Lock (Neo type)

FEATURES

Wrist binder double lock is designed to support, protect and partially immobilize the wrist.

- Extra grip.
- Versatile Velcro closure.
- 4-way compression.

It is made from 3 layered fabric, which is four way stretchable. It provides optimal compression, more comfort, good grip and better durability. It has appealing aesthetics with good color fastness. Wrap around design ensures easy wearing and removal. It allows customized compression, good fitting, better support and wrist stabilization.

Size	Item Code
S	RA5538S
M	RA5538M
L	RA5538L



Elastic Cockup-Splint (Short Type)-Beige

FEATURES

- Pre-shaped mouldable palmar splint allows immobilization in dorsiflexion and flexible dorsal splint limits extension.
- Constructed of high quality coated-fabric with foam pad and sweat absorbing lining for comfort.
- Three sets of buckle and velcro straps allow easy application and provide complete immobilization to the wrist.
- Thumb flap with velcro closure keeps the thumb in slight abduction.
- Designed separately for right and left arms.

INDICATED FOR / APPLICATION

- Strain and sprain in the wrist • Minor fractures of the wrist & distal forearm • Early cast removal



L

Left

Size	Item Code
S	RA5543L.S
M	RA5543L.M
L	RA5543L.L

R

Right

Size	Item Code
S	RA5543R.S
M	RA5543R.M
L	RA5543R.L

Elastic Cockup-Splint (Short Type)-Grey

FEATURES

- Pre-shaped mouldable palmar splint allows immobilization in dorsiflexion and flexible dorsal splint limits extension.
- Constructed of high quality coated-fabric with foam pad and sweat absorbing lining for comfort.
- Three sets of buckle and velcro straps allow easy application and provide complete immobilization to the wrist.
- Thumb flap with velcro closure keeps the thumb in slight abduction.
- Designed separately for right and left arms.

INDICATED FOR / APPLICATION

- Strain and sprain in the wrist • Minor fractures of the wrist & distal forearm • Early cast removal



L

Left

Size	Item Code
S	RA5544L.S
M	RA5544L.M
L	RA5544L.L

R

Right

Size	Item Code
S	RA5544R.S
M	RA5544R.M
L	RA5544R.L

FINGER, WRIST, ELBOW & ARM SPLINTS

Elastic Cockup-Splint (Long Type-Beige)



FEATURES:

For Sprains And Strains In Wrist
Unique Design Provides Firm Support To Wrist
Can Be Used In Left And Right Hand
Easy To Apply And Remove
Comfortable And Quick Pain Relief
Useful In Mild To Moderate Sprains And Strains

L

Left

Size	Item Code
S	RA5547L.S
M	RA5547L.M
L	RA5547L.L

R

Right

Size	Item Code
S	RA5547R.S
M	RA5547R.M
L	RA5547R.L

Elastic Cockup-Splint (Long Type-Beige)



FEATURES:

For Sprains And Strains In Wrist
Unique Design Provides Firm Support To Wrist
Can Be Used In Left And Right Hand
Easy To Apply And Remove
Comfortable And Quick Pain Relief
Useful In Mild To Moderate Sprains And Strains

L

Left

Size	Item Code
S	RA5548L.S
M	RA5548L.M
L	RA5548L.L

R

Right

Size	Item Code
S	RA5548R.S
M	RA5548R.M
L	RA5548R.L

Cockup-Splint with Thumb Support



FEATURES

To immobilize wrist movement, finger & thumb movement also.
Malleable structure can be moulded to fit the contour of for hand and arm

Dorsiflexion angle can be fixed by adjusting velcro straps
Foam padding provides excellent comfort to the patient.

Useful for immobilization in post fracture, tendon injuries of the wrist.

L

Left

Size	Item Code
S	RA5550L.S
M	RA5550L.M
L	RA5550L.L

R

Right

Size	Item Code
S	RA5550R.S
M	RA5550R.M
L	RA5550R.L

Dynamic Cock-up Splint

L

Left

Size	Item Code
S	RA5552L.S
M	RA5552L.M
L	RA5552L.L

R

Right

Size	Item Code
S	RA5552R.S
M	RA5552R.M
L	RA5552R.L

FEATURES

Designed for dynamic splinting following mcp arthroplasty or radial nerve injury or any time extension.

to keep the forearm in desired angular position

Unique spring mechanism allows perfect angle of flexion and extension

Preshaped splint helps to adjust forearm at desired position to achieve quick healing.

Useful in sprains, post fracture, arthritis



Dynamic Radial Palsy Splint with Finger Thumb Attachments

L

Left

Size	Item Code
S	RA5554L.S
M	RA5554L.M
L	RA5554L.L

R

Right

Size	Item Code
S	RA5554R.S
M	RA5554R.M
L	RA5554R.L

FEATURES

- Well padded, unbreakable, washable, reusable and economic.
- Finger-thumb attachments are optional.
- Adjustable wrist-tension by detaching one or two springs.

INDICATED FOR / APPLICATION

- Used for mobilization exercises of fingers and wrist
- For physiotherapy after extensor tendon repair



Volkmann's Turn Buckle Splint

L

Left

Size	Item Code
S	RA5560L.S
M	RA5560L.M
L	RA5560L.L

R

Right

Size	Item Code
S	RA5560R.S
M	RA5560R.M
L	RA5560R.L

FEATURES

To correct contractors of soft tissues

Reverting on plastic material to protect against undue force while correcting

Turn buckle attached to give proper handling of the splints

Very light weight and soft padded for complete comfort

Waterproof, unbreakable design

It is useful in faster correction of contractors of soft tissues



FINGER, WRIST, ELBOW & ARM SPLINTS

Cuff Belt



FEATURES

- Four vertically sewn weight pockets permits cuff to bend freely to conform to any wrist or ankle contour.
- This unique design permits interchangeable use as a wrist or ankle weight cuff.

INDICATED FOR / APPLICATION

- The weight cuff is the only weight belt you will need for a complete progressive exercise or therapy program

Size	Item Code
1/2kg	RA5570A
1kg	RA5570B
2kg	RA5570C
3kg	RA5570D
4kg	RA5570E
5kg	RA5570F

Tennis Elbow Support-Beige



FEATURES

- High strength synthetic strap exerts counter-rotation force at the musculo-tendon region of pronator-supinator muscles (lateral epicondylar region).
- Foam padding offers warmth, compression and support to the inflamed region.
- Velcro with buckle closure enables quick and easy application of the product.

INDICATED FOR / APPLICATION

- Preventive or non-surgical care for Lateral Epicondylitis (Tennis Elbow)

Size	Item Code
Universal	RA5572

Tennis Elbow Support-Grey



FEATURES

- High strength synthetic strap exerts counter-rotation force at the musculo-tendon region of pronator-supinator muscles (lateral epicondylar region).
- Foam padding offers warmth, compression and support to the inflamed region.
- Velcro with buckle closure enables quick and easy application of the product.

INDICATED FOR / APPLICATION

- Preventive or non-surgical care for Lateral Epicondylitis (Tennis Elbow)

Size	Item Code
Universal	RA5573

Tennis Elbow Protector Deluxe, Neo Type

FEATURES

- For relief from tennis elbow, forearm pain.
- Extra foam padding for support to the inflaned region
- High strength strap exerts counter-rotation force at the musculo-tendon region of pronator-supinator.
- Special elastic material ensures strong grip and comfortable compression
- Also recommended in sport injuries.

Size	Item Code
Universal	RA5574

Arm Tourniquet

FEATURES:

- To control or stop blood flow for temperary period
- Easy to apply

Size	Item Code
S(8")	RA5575S
M(10")	RA5575M
L(12")	RA5575L
XL(14")	RA5575XL

Arm Tourniquet

FEATURES:

- To control or stop blood flow for temperary period
- Easy to apply

Size	Item Code
S(8")	RA5576S
M(10")	RA5576M
L(12")	RA5576L
XL(14")	RA5576XL

FINGER, WRIST, ELBOW & ARM SPLINTS

Elastic Tubular Elbow Support Without Centre Hole Delux



FEATURES

Easy to apply a simple pull over application.

Special knitting to cover the affected area in all four ways, which assures effective compression, instant relief, and retains the body heat.

Scientific design ensures uniform gripping, excellent aesthetics and durability.

Special fabric assures better comfort without rash and allergy

Easy to clean, maintain and movement of elbow

It can be used during daily routines, excercises, sports and post surgicals rehabilitation.specially recommended in sprains and strains in elbow joints during sport injuries

It is used to reduce pains, control sprains and strains.provides support to elbow joint and surrounding muscles

Size	Item Code
S	RA5581S
M	RA5581M
L	RA5581L

Elastic Tubular Elbow Support Without Centre Hole premium



FEATURES

Easy to apply. a simple pull over application.

Special knitting to cover the affected area in all four ways, which assures effective compression, instant relief, and retains the body heat.

Scientific design ensures uniform gripping, excellent aesthetics and durability.

Special fabric assures better comfort without rash and allergy

Easy to clean, maintain AND movement of elbow

It can be used during daily routines, excercises, sports and post surgicals rehabilitation.specially recommended in sprains and strains in elbow joints during sport injuries

It is used to reduce pains, control sprains and strains.provides support to elbow joint and surrounding muscles

Size	Item Code
S	RA5583S
M	RA5583M
L	RA5583L

Elastic Tubular Elbow Support Without Centre Hole premium



FEATURES

Easy to apply. a simple pull over application.

Special knitting to cover the affected area in all four ways, which assures effective compression, instant relief, and retains the body heat.

Scientific design ensures uniform gripping, excellent aesthetics and durability.

Special fabric assures better comfort without rash and allergy

Easy to clean, maintain AND movement of elbow

It can be used during daily routines, excercises, sports and post surgicals rehabilitation.specially recommended in sprains and strains in elbow joints during sport injuries

It is used to reduce pains, control sprains and strains.provides support to elbow joint and surrounding muscles

Size	Item Code
S	RA5585S
M	RA5585M
L	RA5585L

Humeral fracture brace

To increase the range of motion of fractured shoulder and elbow
Applied directly on fracture like claw around humerus
Three independent levels of tension control
It is useful in decreasing the rehabilitation time with lessen muscular atrophy

Size	Item Code
S	RA5588S
M	RA5588M
L	RA5588L
XL	RA5588XL



Cuff 'N' Collar

FEATURES

- Made from soft foam in a stockinet sleeve, plastic fastening ties are provided for easy application and removal.

INDICATED FOR / APPLICATION

- Support at wrist and elbow without obscuring the limb under treatment or applying any force to the injured side.

Size	Item Code
Universal	RA 5590



Arm Sling Pouch Adjustable (Deluxe)-Beige

FEATURES

- A high-quality foam fused fabric envelope for the arm provides good, comfortable support and improves patient compliance.
- Non-stretching strap helps to reduce pressure and evenly distribute weight of the affected arm over the other shoulder.
- An adjustable, broad shoulder pad further reduces pressure on neck & shoulder regions.
- Simple velcro closure using plastic buckles.
- Special thumb loop is provided inside the fabric envelope for 'parking' the hand when

Size	Item Code
S	RA5591S
M	RA5591M
L	RA5591L
XL	RA5591XL

INDICATED FOR / APPLICATION

Good support to the arm or shoulder in case of closed or surgical management of strain, sprain and fractures.



FINGER, WRIST, ELBOW & ARM SPLINTS

Arm Sling Pouch Adjustable (Deluxe)-Grey



FEATURES

- A high-quality foam fused fabric envelope for the arm provides good, comfortable support and improves patient compliance.
- Non-stretching strap helps to reduce pressure and evenly distribute weight of the affected arm over the other shoulder.
- An adjustable, broad shoulder pad further reduces pressure on neck & shoulder regions.
- Simple velcro closure using plastic buckles.
- Special thumb loop is provided inside the fabric envelope for 'parking' the hand when desired.

INDICATED FOR / APPLICATION

- Good support to the arm or shoulder in case of closed or surgical management of strain, sprain and fractures.

Size	Item Code
S	RA5592S
M	RA5592M
L	RA5592L
XL	RA5592XL

Arm Sling Pouch Adjustable (Ordinary)



FEATURES

- To provide optimum support to injured, painful arm by reducing neck and shoulder pressure
- Special cotton fabric offers comfortable and soft feeling to injured arm portion
- Easily operative buckle helps in easy application and closer of pouch arm sling
- Thumb loop inside for desired hand rest
- Adjustable special foam pad in shoulder takes load of complete arm
- Very useful in countries with hot and humid weather
- It is useful in fractures and dislocations, sprains and strains, to provide support and immobilize post operative arm

Size	Item Code
S	RA5593S
M	RA5593M
L	RA5593L
XL	RA5593XL

Arm Sling Pouch Adjustable (Premium)-Beige



FEATURES

- Soft quality fabric as envelope for resting the forearm with elbow bent at 90 degree.
- Non-stretching strap helps to reduce pressure & evenly distribute weight of affected arm over the other shoulder.
- Adjustable, broad shoulder pad further reduces pressure on the neck & shoulder region.
- Simple closure, using velcro & plastic buckles.
- Special thumb loop is provided inside the fabric envelope for 'parking' the hand when desired.

INDICATED FOR / APPLICATION

- Support to the arm or shoulder in the case of closed or surgical management of • Strain, sprain • Fractures and dislocations

Size	Item Code
S	RA5594S
M	RA5594M
L	RA5594L
XL	RA5594XL

Arm Sling Pouch Adjustable (Premium)-Grey

FEATURES

- Soft quality fabric as envelope for resting the forearm with elbow bent at 90 degree.
- Non-stretching strap helps to reduce pressure & evenly distribute weight of affected arm over the other shoulder.
- Adjustable, broad shoulder pad further reduces pressure on the neck & shoulder region.
- Simple closure, using velcro & plastic buckles.
- Special thumb loop is provided inside the fabric envelope for 'parking' the hand when desired.


INDICATED FOR / APPLICATION

- Support to the arm or shoulder in the case of closed or surgical management of • Strain, sprain • Fractures and dislocations

Size	Item Code
S	RA5595S
M	RA5595M
L	RA5595L
XL	RA5595XL



KNEE & ANKLE SPLINTS



SIZE CHART

Getting the right splint starts with the right measurement. Our splint are custom made to fit a variety of clients.

- Using a tape measure where indicated arrow.
- Double check the accuracy of the measurement.

* Be sure limb is relaxed during the measurement.



HOW TO MEASURE:

- Measure approx. 6" above to get the thigh circumference.
- Measure 2 to 5" below to get the calf circumference.
- Wrap the measurement tape on ankle to measure the circumference.



Ankle

	XS	S	M	L	XL	XXL	U
Circum. (Inch)	—	7.5-8"	8.5-9.5"	9.5-10.5"	11-12"	11.5-12.5"	Universal
Circum. (cm)	—	19-21	22-24	25-27	28-30	30-32	Universal

Calf

Circum. (Inch)	—	11-13"	12.5-15"	14-16.5"	15.5-18"	16.5-20"	Universal
Circum. (cm)	—	28-34	32-38	36-42	40-46	42-50	Universal

Thigh

Circum. (Inch)	—	16.5-22"	19-25"	21-28"	23.5-30.5"	25.5-33.5"	Universal
Circum. (cm)	—	42-57	48-64	54-71	60-78	65-85	Universal

NOTE: Measurement at 6" above knee center.

- * Customized sizes can be made as per requirement.
- * This size chart is designed to be used as a general guide reference.

Cervical Splints

Shoulder, Rib & Clavicle Splints

Back & Posture Splints

Abdominal Splints

Finger, Wrist, Elbow & Arm Splints

Knee & Ankle Splints

Traction Kits

Orthopaedic General Aids & Gel Pads

Elastic Tubular Knee Support deluxe



FEATURES

Easy to apply, a simple pull over application.

Special knitting to cover the affected area in 2 ways, which assures effective compression, instant relief, and retains the body heat.

Scientific design ensures uniform gripping, excellent aesthetics and durability.

Special fabric assures better comfort.

Easy to clean and maintain.

It can be used during daily routines, exercises, sports and post surgical rehabilitation.

It is used to support knee joint in different conditions like weak knee due to old age, arthritis, sport injuries, support for heavy workouts in sports, post orthopedic surgical rehabilitation

Size	Item Code
S	RA5600S
M	RA5600M
L	RA5600L
XL	RA5600XL

Elastic Tubular Knee Support premium



FEATURES

Easy to apply, a simple pull over application.

Special knitting to cover the affected area in 4 ways, which assures effective compression, instant relief, and retains the body heat.

Scientific design ensures uniform gripping, excellent aesthetics and durability.

Special fabric assures better comfort.

Easy to clean and maintain.

It can be used during daily routines, exercises, sports and post surgical rehabilitation.

It is used to support knee joint in different conditions like weak knee due to old age, arthritis, sport injuries, support for heavy workouts in sports, post orthopedic surgical rehabilitation

Size	Item Code
S	RA5602S
M	RA5602M
L	RA5602L
XL	RA5602XL

Elastic Tubular Knee Support premium Ultra



FEATURES

Easy to apply, a simple pull over application.

Special knitting to cover the affected area in 4 ways, which assures effective compression, instant relief, and retains the body heat.

100 % NYLON outside and COTTON inside.

Scientific design ensures uniform gripping, excellent aesthetics and durability.

Special fabric assures better comfort.

Easy to clean and maintain.

It can be used during daily routines, exercises, sports and post surgical rehabilitation.

It is used to support knee joint in different conditions like weak knee due to old age, arthritis, sport injuries, support for heavy workouts in sports, post orthopedic surgical rehabilitation

Size	Item Code
S	RA5604S
M	RA5604M
L	RA5604L
XL	RA5604XL

Elastic Tubular Knee Support With Centre Hole-delux

FEATURES

Easy to apply, a simple pull over application.
Special knitting to cover the affected area in 2 ways, which assures effective compression, instant relief, and retains the body heat.
Open Patella hole.
Scientific design ensures uniform gripping, excellent aesthetics and durability.
Special fabric assures better comfort.
Easy to clean and maintain.
It can be used during daily routines, exercises, sports and post surgical rehabilitation.

It is used to support knee joint in different conditions like weak knee due to old age, arthritis, sport injuries, support for heavy workouts in sports, post orthopedic surgical rehabilitation

Size	Item Code
S	RA5606S
M	RA5606M
L	RA5606L
XL	RA5606XL



Elastic Tubular Knee Support With Centre Hole premium

FEATURES

Easy to apply, a simple pull over application.
Special knitting to cover the affected area in 4 ways, which assures effective compression, instant relief, and retains the body heat.
Scientific design ensures uniform gripping, excellent aesthetics and durability.
Open Patella hole.
Special fabric assures better comfort.
Easy to clean and maintain.
It can be used during daily routines, exercises, sports and post surgical rehabilitation.

It is used to support knee joint in different conditions like weak knee due to old age, arthritis, sport injuries, support for heavy workouts in sports, post orthopedic surgical rehabilitation

Size	Item Code
S	RA5608S
M	RA5608M
L	RA5608L
XL	RA5608XL



Elastic Tubular Varicose Stockings-below Knee-deluxe

FEATURES

Easy to apply, a simple pull over application.
Special knitting to cover the affected area in 2 ways, which assures effective compression, instant relief, and retains the body heat.
To relieve discomfort from tired swollen legs
Provides support to thromboembolism
Seamless stocking from below knee to ankle
Scientific design ensures uniform gripping, excellent aesthetics and durability.
Special fabric assures better comfort.
Easy to clean and maintain.

It can be used during daily routines, exercises, sports and post surgical rehabilitation.

Size	Item Code
S	RA5610S
M	RA5610M
L	RA5610L
XL	RA5610XL



Elastic Tubular Varicose Stockings-thigh Length-deluxe



FEATURES

Easy to apply, a simple pull over application.
 Special knitting to cover the affected area in 2 ways, which assures effective compression, instant relief, and retains the body heat.
 To relieve discomfort from tired swollen legs
 Provides support to thromboembolism
 Seamless stocking from mid thigh to ankle
 Scientific design ensures uniform gripping, excellent aesthetics and durability.
 Special fabric assures better comfort.
 Easy to clean and maintain.
 It can be used during daily routines, exercises, sports and post surgical rehabilitation.

Size	Item Code
S	RA5612S
M	RA5612M
L	RA5612L
XL	RA5612XL

Elastic Tubular Varicose Stockings-below Knee-premium



FEATURES

Easy to apply, a simple pull over application.
 Special knitting to cover the affected area in 4 ways, which assures effective compression, instant relief, and retains the body heat.
 To relieve discomfort from tired swollen legs
 Provides support to thromboembolism
 Seamless stocking from below knee to ankle
 Scientific design ensures uniform gripping, excellent aesthetics and durability.
 Special fabric assures better comfort.
 Easy to clean and maintain.
 It can be used during daily routines, exercises, sports and post surgical rehabilitation.

Size	Item Code
S	RA5614S
M	RA5614M
L	RA5614L
XL	RA5614XL

Elastic Tubular Varicose Stockings-thigh Length-premium



FEATURES

Easy to apply, a simple pull over application.
 Special knitting to cover the affected area in 4 ways, which assures effective compression, instant relief, and retains the body heat.
 To relieve discomfort from tired swollen legs
 Provides support to thromboembolism
 Seamless stocking from mid thigh to ankle
 Scientific design ensures uniform gripping, excellent aesthetics and durability.
 Special fabric assures better comfort.
 Easy to clean and maintain.
 It can be used during daily routines, exercises, sports and post surgical rehabilitation.

Size	Item Code
S	RA5616S
M	RA5616M
L	RA5616L
XL	RA5616XL

Elastic Tubular Knee Support With Open Patella Gel-premium

FEATURES

Easy to apply, a simple pull over application.
Special knitting to cover the affected area in 4 ways, which assures effective compression, instant relief, and retains the body heat.
Scientific design ensures uniform gripping, excellent aesthetics and durability.
Open Patella hole with Gel support for more comfort and support at patella.
Special fabric assures better comfort.
Easy to clean and maintain.
It can be used during daily routines, exercises, sports and post surgical rehabilitation.

It is used to support knee joint in different conditions like weak knee due to old age, arthritis, sport injuries, support for heavy workouts in sports, post orthopedic surgical rehabilitation

Size	Item Code
S	RA5620S
M	RA5620M
L	RA5620L
XL	RA5620XL



Knee Support Open Patella Neo With Hinges

FEATURES

It is used in painful knee movement of arthritis & to prevent mild strain & sprains that may occur due to activities of daily living. Excellent quality three layered drytex fabric of 12" width forms the main body of the support.

Anterior patella hole & posterial cut out are provided to avoid bunching

Unique metal splints provide lateral support while controlling movement of the knee

Product design has Velcro closure to fit support around knee.

Size	Item Code
S	RA5622S
M	RA5622M
L	RA5622L
XL	RA5622XL



Knee Support Open Patella Neo Without Hinges

FEATURES

It is used in painful knee movement of arthritis & to prevent mild strain & sprains that may occur due to activities of daily living. Excellent quality three layered drytex fabric of 12" width forms the main body of the support.

Anterior patella hole & posterial cut out are provided to avoid bunching

Product design has Velcro closure to fit support around knee.

Size	Item Code
S	RA5624S
M	RA5624M
L	RA5624L
XL	RA5624XL



KNEE & ANKLE SPLINTS

Rom Brace



FEATURES

Provides in the fracture of knee pre fabricated.
To control and restrict flexion-extension of the knee joint
Scientific design allows controlled movements in any two angles
Can be adjusted easily by changing the position of pin
Knee joint can be immobilized at any angle
Marked angles at metal disc with pin locking for added safety
It can be used for maximum orthopaedic problems related to knee joints.

Size	Item Code
S	RA5626S
M	RA5626M
L	RA5626L
XL	RA5626XL

Elastic Knee Support-ordinary



FEATURES

It is used in painful knee movement of arthritis & to prevent mild strain & sprains that may occur due to activities of daily living.
Fine quality elastic of 9" width forms the main body of the support.
Product design has Velcro closure to fit support around knee.

Size	Item Code
S	RA5630S
M	RA5630M
L	RA5630L
XL	RA5630XL

Elastic Knee Support-delux



FEATURES

It is used in painful knee movement of arthritis & to prevent mild strain & sprains that may occur due to activities of daily living.
Excellent quality elastic of 9" width forms the main body of the support.
Product design has Velcro closure to fit support around knee.

Size	Item Code
S	RA5632S
M	RA5632M
L	RA5632L
XL	RA5632XL

Elastic Knee Support With Hinges Delux

FEATURES

It is used in painful knee movement of arthritis & to prevent mild strain & sprains that may occur due to activities of daily living. Excellent quality elastic of 9" width forms the main body of the support.

Anterior patella hole & posterior cut out are provided to avoid bunching

Unique metal splints provide lateral support while controlling movement of the knee

Product design has Velcro closure to fit support around knee.

Size	Item Code
S	RA5634S
M	RA5634M
L	RA5634L
XL	RA5634XL



Elastic Tubular Knee Support With Hinges(single)-deluxe

FEATURES

Easy to apply, a simple pull over application.

Special knitting to cover the affected area in 2 ways, which assures effective compression, instant relief, and retains the body heat. Unique flexible splints provide lateral support while controlling movement of the knee

Product design has Velcro closure to fit support around knee.

Scientific design ensures uniform gripping, excellent aesthetics and durability.

Special fabric assures better comfort.

Easy to clean and maintain.

It can be used during daily routines, exercises, sports and post surgical rehabilitation.

It is used to support knee joint in different conditions like weak knee due to old age, arthritis, sport injuries, support for heavy workouts in sports, post orthopedic surgical rehabilitation

Size	Item Code
S	RA5636S
M	RA5636M
L	RA5636L
XL	RA5636XL



Elastic Tubular Knee Support With Hinges(single)-premium

FEATURES

Easy to apply, a simple pull over application.

Special knitting to cover the affected area in 4 ways, which assures effective compression, instant relief, and retains the body heat. Unique flexible splints provide lateral support while controlling movement of the knee

Product design has Velcro closure to fit support around knee.

Scientific design ensures uniform gripping, excellent aesthetics and durability.

Special fabric assures better comfort.

Easy to clean and maintain.

It can be used during daily routines, exercises, sports and post surgical rehabilitation.

It is used to support knee joint in different conditions like weak knee due to old age, arthritis, sport injuries, support for heavy workouts in sports, post orthopedic surgical rehabilitation

Size	Item Code
S	RA5638S
M	RA5638M
L	RA5638L
XL	RA5638XL



Knee Brace Short Type-14" Beige



FEATURES

It is used in fractures, ligament or muscle related injuries and dislocation of knee joint
Easy to apply.
Light weight and comfortable.
Short type- 14" in length with different sizes
Covers maximum portion, knee which assures immobilization of affected knee joint.
Scientific design ensures uniform gripping, excellent aesthetics and durability with rigid malleable aluminium splinting at sides.
Special fabric assures better comfort without rash and allergy
Easy to clean and maintain
It can be used during daily routines, post surgical rehabilitation.
Specially recommended in serious knee problems.

Size	Item Code
S	RA5640S
M	RA5640M
L	RA5640L
XL	RA5640XL

Knee Brace Short Type-14" Grey



FEATURES

It is used in fractures, ligament or muscle related injuries and dislocation of knee joint
Easy to apply.
Light weight and comfortable.
Short type- 14" in length with different sizes
Covers maximum portion, knee which assures immobilization of affected knee joint.
Scientific design ensures uniform gripping, excellent aesthetics and durability with rigid malleable aluminium splinting at sides.
Special fabric assures better comfort without rash and allergy
Easy to clean and maintain
It can be used during daily routines, post surgical rehabilitation.
Specially recommended in serious knee problems.

Size	Item Code
S	RA5641S
M	RA5641M
L	RA5641L
XL	RA5641XL

Knee Brace Long Type-19" Beige



FEATURES

It is used in fractures, ligament or muscle related injuries and dislocation of knee joint
Light weight and comfortable.
Long type- 19" in length with different sizes
Covers maximum portion, knee and calf which assure immobilization of affected knee joint.
Scientific design ensures uniform gripping, excellent aesthetics and durability with rigid malleable aluminium splinting at sides.
The stays are removable and moldable for proper fitting.
Buckle & velcro closures above & below knee provide good fit.
It can be used during daily routines, post surgical rehabilitation.
Specially recommended in serious knee problems.

Size	Item Code
S	RA5642S
M	RA5642M
L	RA5642L
XL	RA5642XL

KNEE & ANKLE SPLINTS

Knee Brace Long Type-19" Grey

FEATURES

It is used in fractures, ligament or muscle related injuries and dislocation of knee joint
Light weight and comfortable.
Long type- 19" in length with different sizes
Covers maximum portion, knee and calf which assure immobilization of affected knee joint.
Scientific design ensures uniform gripping, excellent aesthetics and durability with rigid malleable aluminium splinting at sides.
The stays are removable and moldable for proper fitting.
Buckle & velcro closures above & below knee provide good fit.
It can be used during daily routines, post surgical rehabilitation.
Specially recommended in serious knee problems.

Size	Item Code
S	RA5643S
M	RA5643M
L	RA5643L
XL	RA5643XL



Knee Brace Deluxe / Knee Immobilizer-23" Beige

FEATURES

It is used in fractures, ligament or muscle related injuries and dislocation of knee joint
Light weight and comfortable.
Knee immobilizer 23" in length with different sizes for complete immobilization of Knee.
Covers maximum portion, knee and calf which assure immobilization of affected knee joint.
Scientific design ensures uniform gripping, excellent aesthetics and durability with rigid malleable aluminium splinting at sides.
The stays are removable and moldable for proper fitting.
Buckle & velcro closures above & below knee provide good fit.
It can be used during daily routines, post surgical rehabilitation.
Specially recommended in serious knee problems.

Size	Item Code
S	RA5644S
M	RA5644M
L	RA5644L
XL	RA5644XL



Knee Brace Deluxe / Knee Immobilizer-23" Grey

FEATURES

It is used in fractures, ligament or muscle related injuries and dislocation of knee joint
Light weight and comfortable.
Knee immobilizer 23" in length with different sizes for complete immobilization of Knee.
Covers maximum portion, knee and calf which assure immobilization of affected knee joint.
Scientific design ensures uniform gripping, excellent aesthetics and durability with rigid malleable aluminium splinting at sides.
The stays are removable and moldable for proper fitting.
Buckle & velcro closures above & below knee provide good fit.
It can be used during daily routines, post surgical rehabilitation.
Specially recommended in serious knee problems.

Size	Item Code
S	RA5645S
M	RA5645M
L	RA5645L
XL	RA5645XL



KNEE & ANKLE SPLINTS

Thigh Band



FEATURES

- Support in case of quadriceps muscle strain & hamstring pull.
- High quality 9" elastic band provides compression & supports to thigh muscles
- High tenacity strap along the length of product maintains its tubular structure
- Velcro closure easy to use.

Size	Item Code
S	RA5650S
M	RA5650M
L	RA5650L
XL	RA5650XL

Ankle Support



FEATURES

- Can be fixed firmly to the ankle by Velcro straps.
- Constructed of semi-rigid outer shell and soft foam padding.
- Provides lateral support limiting inversion and eversion while allowing normal flexion.
- Designed to be worn comfortably inside laced shoes.
- Velcro closure allows quick application and proper fit.

INDICATED FOR / APPLICATION

- Post-operative support • Sprains/strains • Fractures and cast removal

Size	Item Code
Universal	RA5660

Ankle Brace



FEATURES

- Designed to be worn comfortably inside laced shoe.
- Controls abduction and adduction while allowing flexion and extension.
- Provides efficient support without stays.
- Lace-up design facilitates custom fit.

INDICATED FOR / APPLICATION

- Prevention of ankle injuries during sports and normal activities • To protect ankle in the case of ligament tears, and hairline fractures • Post-operative use and early cast removal

Size	Item Code
S	RA5662S
M	RA5662M
L	RA5662L
XL	RA5662XL

Ankle Brace

FEATURES

It is used in ligament tears & hairline fracture & early cast removal
Durable quoted fabric with inner lining of form comfort & fit around the ankle joint
High quality semi rigid plastic stays on sides of base control abduction.
Restrict flexion-extension movements
Special Velcro straps and plastic buckles assure effective reverse locking with secured fitting to ankle joint.
It can be worn on either left or right ankle

Size	Item Code
S	RA5663S
M	RA5663M
L	RA5663L
XL	RA5663XL



Elastic Ankle Binder-Beige

FEATURES

It is used in a mild ankle sprain during the activities of daily living.
Special elastic fabric in figure eight wrap assures controlled gripping all around the ankle joint and pain relief.
Convenient to apply and can be used for prolonged period
It is used in different conditions like weak ankle due to old age, serious ankle injuries, post operative orthopedic surgical rehabilitation.
Can be worn inside shoes.

Size	Item Code
S	RA5666S
M	RA5666M
L	RA5666L
XL	RA5666XL



Elastic Ankle Binder-Black

FEATURES

It is used in a mild ankle sprain during the activities of daily living.
Special elastic fabric in figure eight wrap assures controlled gripping all around the ankle joint and pain relief.
Convenient to apply and can be used for prolonged period
It is used in different conditions like weak ankle due to old age, serious ankle injuries, post operative orthopedic surgical rehabilitation.
Can be worn inside shoes.

Size	Item Code
S	RA5668S
M	RA5668M
L	RA5668L
XL	RA5668XL



KNEE & ANKLE SPLINTS

Elastic Tubular Anklet With Binder-Grey

FEATURES

It is composed of two components compression ankle sleeve along with ankle binder

Four way stretchable fabric, ensures effective compression and comfort

Special elastic fabric in figure eight wrap assures controlled gripping all around the ankle joint and pain relief.

Convenient to apply and can be used for prolonged period

It is used in different conditions like weak ankle due to old age, serious ankle injuries, and post operative orthopedic surgical rehabilitation.

It is used in a mild ankle sprain during the activities of daily living.

Size	Item Code
S	RA5669S
M	RA5669M
L	RA5669L
XL	RA5669XL



Elastic Tubular Anklet With Binder-Beige

FEATURES

It is composed of two components compression ankle sleeve along with ankle binder

Four way stretchable fabric, ensures effective compression and comfort

Special elastic fabric in figure eight wrap assures controlled gripping all around the ankle joint and pain relief.

Convenient to apply and can be used for prolonged period

It is used in different conditions like weak ankle due to old age, serious ankle injuries, and post operative orthopedic surgical rehabilitation.

It is used in a mild ankle sprain during the activities of daily living.

Size	Item Code
S	RA5670S
M	RA5670M
L	RA5670L
XL	RA5670XL



Elastic Tubular Anklet-deluxe

FEATURES

Easy to apply, a simple pull over application.

It is used to support Ankle joints in different conditions.

Special knitting to cover the affected area in 2 ways, which assures effective compression, instant relief, and retains the body heat.

Scientific design ensures uniform gripping, excellent aesthetics and durability.

Special fabric assures better comfort.

Easy to clean and maintain.

It can be used during daily routines, exercises, sports and post surgical rehabilitation.

It is used to support knee joint in different conditions like weak knee due to old age, arthritis, sport injuries, support for heavy workouts in sports, post orthopedic surgical rehabilitation

Size	Item Code
S	RA5672S
M	RA5672M
L	RA5672L
XL	RA5672XL



KNEE & ANKLE SPLINTS

Elastic Tubular Anklet-premium



FEATURES

Easy to apply, a simple pull over application.

It is used to support Ankle joints in different conditions.

Special knitting to cover the affected area in 4 ways, which assures effective compression, instant relief, and retains the body heat.

Scientific design ensures uniform gripping, excellent aesthetics and durability.

Special fabric assures better comfort.

Easy to clean and maintain.

It can be used during daily routines, exercises, sports and post surgical rehabilitation.

It is used to support knee joint in different conditions like weak knee due to old age, arthritis, sport injuries, support for heavy workouts in sports, post orthopedic surgical rehabilitation

Size	Item Code
S	RA5674S
M	RA5674M
L	RA5674L
XL	RA5674XL

Foot Drop Splint-leaf Type



features

For stabilization of ankle and foot in foot drop condition

Pre-molded, light weight and thermoplastic design

Light weight padded splint to give comfortable wear

Can wear inside shoes.

Velcro straps provided to give proper fitting

It is useful while walking.

Size	Item Code
S	RA5680S
M	RA5680M
L	RA5680L

Foot Drop Splint



FEATURES

- Foot portion ends below the metatarsal heads.
- Can be used as a night splint for foot drop.

INDICATED FOR / APPLICATION

- Lightweight, padded posterior leg splints for temporary or short-term immobilization of fractures and sprains.

Size	Item Code
S	RA5682S
M	RA5682M
L	RA5682L
XL	RA5682XL

Foot / Ankle Leg Splint

L

Left

Size	Item Code
S	RA5684L.S
M	RA5684L.M
L	RA5684L.L
XL	RA5684L.XL

FEATURES

- It is well padded, unbreakable, washable, reusable and economic.

INDICATED FOR / APPLICATION

- It is the most ideal alternate to wire splint or cumbersome POP slab, and is very easy to wear • It allows dressing to be done without disturbing the splint

R

Right

Size	Item Code
S	RA5684R.S
M	RA5684R.M
L	RA5684R.L
XL	RA5684R.XL

L

Left

Size	Item Code
1	RA5687L.1
2	RA5687L.2
3	RA5687L.3
4	RA5687L.4
5	RA5687L.5
6	RA5687L.6
7	RA5687L.7

FEATURES

- Light weighted and comfortable.
- Weight bearing is allowable depending upon the condition of bone/fracture.
- Cosmetically highly accepted.
- Removable and hence care of skin is easy.

INDICATED FOR / APPLICATION

- Period in plaster cast can be reduced significantly.

R

Right

Size	Item Code
1	RA5687R.1
2	RA5687R.2
3	RA5687R.3
4	RA5687R.4
5	RA5687R.5
6	RA5687R.6
7	RA5687R.7

Tibial Functional Brace

Cuff Belt

FEATURES


- Four vertically sewn weight pockets permits cuff to bend freely to conform to any wrist or ankle contour.
- This unique design permits interchangeable use as a wrist or ankle weight cuff.

INDICATED FOR / APPLICATION

- The weight cuff is the only weight belt you will need for a complete progressive exercise or therapy program.

Size	Item Code
1/2kg	RA5570A
1kg	RA5570B
2kg	RA5570C
3kg	RA5570D
4kg	RA5570E
5kg	RA5570F

TRACTION KITS



Size Chart

Cervical Splints

Shoulder, Rib & Clavicle Splints

Back & Posture Splints


Abdominal Splints

Finger, Wrist, Elbow & Arm Splints

Knee & Ankle Splints

Traction Kits

Orthopaedic General Aids & Gel Pads



TRACTION KITS

Pelvic Traction Kit (Complete Set)



FEATURES

- Foam padded, pelvic traction belt is designed to supply traction to the lower spine, including the lumbar and sacroiliac areas.

INDICATED FOR / APPLICATION

- Gives traction of the lumbosacral region for relief of low back pain and to maintain normal alignment of the vertebrae.

Kit Content

- Pelvic Belt
- Traction Pulley
- Heavy Spread Bar,
- String
- Sand or Water Bag

Size	Item Code
S	RA5700S
M	RA5700M
L	RA5700L
XL	RA5700XL

Leg Traction Kit (Complete Set)



FEATURES

- To support the traction arrangement for lower back pain, knee region
- Special nonslip foam lined stirrup makes it useful for sensitive skin
- With heavy duty straps and strong metal rings belt is directly connectable to traction kit, Traction Pulley & Sand Bag
- Easy to apply and remove
- Gives controlled and uniform traction
- It is useful in lower back and knee related pains

Kit Content

- Leg Traction Belt
- Traction Pulley
- Spread Bar
- String
- Sand or Water Bag

Size	Item Code
S	RA5710S
M	RA5710M
L	RA5710L
XL	RA5710XL

Foot/Skin Traction Kit (Complete Set)

FEATURES

- Foot/skin traction kit is used as a below knee traction.
- The skin traction strips provide maximum aeration to the skin while securely attaching its center extremity to the spreader.

Kit Content

- 3 Layered Foam Fabri
- M-Angle,
- String
- Cotton Tape



Size	Item Code
L	RA5720L

Ankle Traction Kit (Complete Set)

Features

- To help in traction of ankle lower leg section with heavy duty straps and strong rings belt is directly connectable to traction kit, Traction Pulley & Sand Bag
- Easy to apply and remove
- Gives controlled and uniform traction.
- It is useful for traction in ankle region

Kit Content

- Ankle Belt
- Traction Pulley
- String
- Sand or Water Bag



Size	Item Code
S	RA5730S
M	RA5730M
L	RA5730L
XL	RA5730XL

TRACTION KITS

Cervical Traction Kit (Sleeping Type) Complete Set



FEATURES

- This is portable apparatus designed for use at home, private clinics and hospitals.
- Very useful where continuous traction is recommended by doctors.

INDICATED FOR

- Sleeping traction is used to treat cervical spine

Kit Content

- Traction Header
- Traction Pulley
- Spread Bar
- String
- Sand or Water Bag

RA 5740

Available Size

One size fits all adults

Cervical Traction Kit (Sitting Type) Complete Set



FEATURES

- This is portable apparatus designed for use at home, private clinics and hospitals.
- Longitudinal stretching by traction effectively relieves pain associated with cervical spondylosis.

INDICATED FOR / APPLICATION

- Cervical spine traction in sitting position to relieve pain associated with muscle spasm (contraction due to tension, surgery etc.)

Kit Content




- Traction Header
- Door Traction Pulley
- Spread Bar
- String
- Sand or Water Bag


RA 5750

Available Size

One size fits all adults

GENERAL AIDS





Cervical Splints

Shoulder, Rib & Clavicle Splints

Back & Posture Splints

Abdominal Splints

Finger, Wrist, Elbow & Arm Splints

Knee & Ankle Splints

Traction Kits

Orthopaedic General Aids & Gel Pads

Orthopaedic Back Rest Regular



FEATURES

Orthopaedically designed for sitting comfort.

It gives spinal support when placed on chair, auto seat & wheel chair.

It is ideal for slip disc patients.reduces discomfort and pain after surgery.

It improves posture when placed on chairs.

to maintain the natural curve of lower back and

Helps in releasing the pressure on vertibras and gives correct body posture while sitting

Size	Item Code
Universal	RA 5800

Orthopaedic Back Rest short type



FEATURES

- Lightweight and portable back rest which gives spinal support when placed on any chair, car seat, sofa, wheel chair, bed or on the floor.
- This is a brace, which you do not have to wear.
- It doesn't cure back ache but eases the pain associated with long hours of chair sitting.

INDICATED FOR / APPLICATION

- It is ideal for slip disc patients • Removes discomfort and pain after surgery

Size	Item Code
Universal	RA 5802

Orthopaedic Back Rest Long Type

FEATURES

Orthopaedically designed for sitting comfort.
It gives spinal support when placed on chair, auto seat & wheel chair.
It is ideal for slip disc patients.reduces discomfort and pain after surgery.
It improves posture when placed on chairs.
To maintain the natural curve of lower back and
Helps in releasing the pressure on vertibras and gives correct body posture while sitting



Size	Item Code
Universal	RA 5804

Orthopaedic Back Rest Long Type with neck Support

FEATURES

- Lightweight and portable back rest which gives spinal support when placed on any chair, car seat, sofa, wheel chair, bed or on the floor.
- This is a brace, which you do not have to wear.
- It doesn't cure back ache but eases the pain associated with long hours of chair sitting.

INDICATED FOR / APPLICATION

- It is ideal for slip disc patients • Removes discomfort and pain after surgery



Size	Item Code
Universal	RA 5806

Orthopaedic Coccyx Wedge Seat



FEATURES:

Scientific design release pressure points on coccyx region
Special foam provides total comfort to complete lower back section
It is useful to prevent severe back pain in prolonged sitting position.
To release the pressure on lower spinal discs and tail bone(coccyx) due to prolonged sitting



Size	Item Code
Universal	RA 5810

Orthopaedic Coccyx Wedge Seat-Ring Cushion



FEATURES:

This air seat cushion is especially designed for wheel chair and daily prolonged sitting use.

- Durable material and heat sealed construction to prevents air from leaking.
- Light weight and foldable for carrying around with you.
- Helps to relieve pain and stress when sitting.
- Air inflatable waffle seat cushion heat seated construction for durability.



Size	Item Code
Universal	RA 5811

GENERAL AIDS

Orthopaedic Heating Pad

FEATURES

To provide relief from pain for arthritis patient, decrease joint stiffness

Three different levels of heating controlled by sliding switch to produce appropriate heat as per patients requirement
Special outer cotton cover keeps the belt in touch with patient and transfers the maximum heat to the patient.

Use: in back, thigh, shoulder pains



Size	Item Code
Universal	RA5820

Cast Shoes

FEATURES

- An ideal walking aid for patients, whose leg is in plaster. Prevents the plaster from chipping, cracking or soiling.
- Facilitates stability while walking.



Size	Item Code
S	RA5830S
M	RA5830M
L	RA5830L
XL	RA5830XL

Exercise Ball Neuro & Ortho (Gel) - Soft



FEATURES

- Gel material.
- To increase the strength of intrinsic muscles of hand.
- It is helpful in fingers & hands exercise for grip strength.
- It gives mobility & fine motor skill.

Size	Item Code
Universal	RA 5840

Exercise Ball Neuro & Ortho (Gel) - Hard



FEATURES

- Gel material.
- To increase the strength of intrinsic muscles of hand.
- It is helpful in fingers & hands exercise for grip strength.
- It gives mobility & fine motor skill.

Size	Item Code
Universal	RA 5841

Exercise Ball Neuro & Ortho (Foam)



FEATURES

- Soft material.
- To increase the strength of intrinsic muscles of hand.
- It is helpful in fingers & hands exercise for grip strength.
- It gives mobility & fine motor skill.

Size	Item Code
Universal	RA 5842

GENERAL AIDS

Finger Exerciser Plastic

FEATURES

- For Hand exercises.
- To increase the strength of intrinsic muscles of hand.
- Can effectively the finger blood circulation enhanced finger strength, particularly effective.
- Light weight, pre-molded with sliding frame work.

Size	Item Code
Universal	RA5850



RUBBER Cast HEEL

FEATURES

Used with weight bearing foot cast
Moulded out of light weight durable & resilient polymer
Durable rubber with graduations to provides resistance against slipping
Arch provided to give proper stability



Size	Item Code
S	RA5860S
M	RA5860M
L	RA5860L

C.T.E.V SPLINT

FEATURES

For correcting mild ctev deformities
Features and benefits:
Washable, light weight design
Smooth padding provides extra comfort to infants
Useful to avoid repeated plastering and surgery in infants for ctev deformities



Size	Item Code
1	RA5870A
2	RA5870B
3	RA5870C
4	RA5870D
5	RA5870E
6	RA5870F
7	RA5870G



Cotton Crepe Bandage B.P with Fast Edges

FEATURES:

In the treatment of sprains and dislocated painful joints.
Excellent elasticity provides uniform gripping and compression
Easy to apply and remove.
Multiple use in problems related to sprains and strains

Size	Item Code
06x4cm Stretched	RA5880A
08x4cm Stretched	RA5880B
10x4cm Stretched	RA5880C
15x4cm Stretched	RA5880D

Orthopaedic Stockinet Roll (RIB Cloth)

FEATURES

Inside the plaster
Made from very soft non absorbent tubular rib weaved cotton
Easy to apply

Size	Item Code
2.5"x20mts	RA5890A
3.0"x20mts	RA5890B
4.0"x20mts	RA5890C
5.0"x20mts	RA5890D
6.0"x20mts	RA5890E

Thomas Splint

FEATURES

Leather / rexine covered swivel ring, tinned, used for fractures of upper & lower thigh.


Size	Item Code
XS	RA5900XS
S	RA5900S
M	RA5900M
L	RA5900L
XL	RA5900XL

Cramer's Wire Splint

FEATURES

It is used as emergency splint before applying the cast
Made of moldable wire and vinyl coated for a soft pliable surface finish
Light weight and malleable

Available



Size		Item Code
2"	X 24"	RA5910A
3"	X 36"	RA5910B
4"	X 36"	RA5910C
5"	X 36"	RA5910D
6"	X 36"	RA5910E

GENERAL AIDS

GEL PADS

Toe Separator



FEATURES

To relieve friction between overlapping and pinch toes
Scientific design helps in spreading and aligning the toes along with absorption of pressure.
Easy to use inside shoes
It is useful in pain and soreness

Size	Item Code
Set of 2	RA 6000
Set Of 4	RA 6001

Arch Pad (Gel)



FEATURES

To provide additional arch support to maintain proper foot position.
Special gel supports fallen arches weak arches high arches.
Useful in common orthopaedic problems of arch region

Size	Item Code
S	RA6011S
M	RA6011M
L	RA6011L

Heel Cushion (Gel)



FEATURES

To provide cushioning effect to heel against the painful un-operatable over growth and development of calcaneal spur
Heel pad made of gel
Scientific design provides comfortable cushioning and adjusts automatically as per the pressure on toe region which gives relief from pain to the user
Durable and long life
Easy to clean and maintain
It is useful in common orthopaedic problems of heel region

Size	Item Code
M	RA6021M
L	RA6021L

GENERAL AIDS GEL PADS

Heel Pad Female (Gel)

PRODUCT FEATURES:

To provide cushioning effect to heel against the painful un-operatable over growth and development of calcaneal spur
Heel pad female made of gel
Scientific design provides comfortable cushioning and adjusts automatically as per the pressure on toe region which gives relief from pain to the user
Durable and long life
Easy to clean and maintain
It is useful in common orthopaedic problems of heel region

Size	Item Code
Standard	RA 6031



Heel Pad (Gel)

FEATURES

To provide cushioning effect to heel against the painful un-operatable over growth and development of calcaneal spur
Heel pad made of gel
Scientific design provides comfortable cushioning and adjusts automatically as per the pressure on toe region which gives relief from pain to the user
Durable and long life
Easy to clean and maintain
It is useful in common orthopaedic problems of heel region

Size	Item Code
Standard	RA 6035



Heel Pad (Donut)

FEATURES

To provide cushioning effect to heel against the painful un-operatable over growth and development of calcaneal spur
Scientific design provides comfortable cushioning and adjusts automatically as per the pressure on toe region which gives relief from pain to the user
Durable and long life
Easy to clean and maintain
It is useful in common orthopaedic problems of heel region



Size	Item Code
Standard	RA 6040

GENERAL AIDS GEL PADS

Full Length Insole (Gel)



FEATURES

Made of gel

For complete heel to toe protection and comfort

Absorbs maximum shock impact on heel to toe region

Effective in foot as well as back pain

Easy to use

Long life

Size	Item Code
S	RA6045S
M	RA6045M
L	RA6045L

TERMS & CONDITIONS

- All sizes, dimensions, values etc. in this catalogue are approximate. We reserve the right to upgrade / improve / alter the specifications of our products without notice.
- We reserve the right to supply different colour/design than those shown in the catalogue, which depends on the availability and trend of raw material at the time of supply
- Use of all our products is restricted to professionals who are qualified and trained to use them. We are not responsible for any damage which may occur due to improper or unauthorized use. In case of malfunctioning of the product, the company's liability is limited to the maximum extent of the cost of equipment.
- It is users' responsibility to keep the products and packing away from the reach of unauthorised people, children and animals and dispose of them in an environment friendly manner.
- You must keep proper record of your sales. This will assist product re-calls (if required).

Rehabilitation Aids & Orthopaedic Soft Goods



Henley on Thames

RG91NN

United Kingdom

Tel: +44 7788801854

email: info@orthotraumaint.co.uk

www.orthotraumaint.co.uk

