

Face Mask


Disposable Protective Face Mask
3Ply



- Spandex Loop Non Latex
- Nose Pin for Comfort
- Made by Automatic Plant
- Automatic Welded Ear Loop
- Available in different options of filter protection & fluid resistance
- Fine quality fabric comfort to inhale and exhale
- We maintain hygiene in manufacturing process
- Test results by Certified Testing Institute available on request
- Class I CE Certificate of conformity & ISO 9001:2015

Salient Features

3 Layers of Protection



Nose Pin
for Comfort



Made By
Automatic Plant



Different options
of **Filter Protection**
& **Fluid Resistance**

Level 1 Mask	Level 2 Mask	Level 3 Mask
<p>Item Code: HC406.3N & HC406.3NB</p>	<p>Item Code: HC406.3S & HC406.3SB</p>	<p>Item Code: HC406.3F & HC406.3FB</p>
<ul style="list-style-type: none"> • Designed for moderate risk fluid environments. • Provides medium bacterial and particle filtration. • Great comfort and breathability. • Made of 3 layers of PP Spunbond Nonwoven material. 	<ul style="list-style-type: none"> • Designed for moderate to higher risk fluid environments. • Provides good bacterial and particle filtration. • Good balance of performance and comfort. • Made of 2 layers of PP Spunbond Nonwoven & 1 layer of SMS material. 	<ul style="list-style-type: none"> • Designed for higher risk fluid environments. • Provides excellent bacterial and particle filtration. • Maintain breathability while providing a high level of fluid protection. • Made of 2 layers of PP Spunbond Nonwoven & 1 layer of Melt Blown material.