

N-95

Protective Mask



5 Layers
Protection

Above **95%**
Filtration Efficiency



High
Protection



Comfortable
& Breathable



Anti-Dust
Anti-Odor



Infection
Prevention

N-95 Protective Mask

FEATURES



Adjustable
Nose Bridge




Ultrasonic
Stitching



5 Layers
Protection



Elastic
Ear Loops



3D Middle Half
Fold Design

USAGE

Wearing a face mask will protect the individual from contracting any airborne infectious germs and prevent the spread of infection. When someone coughs, talks & sneezes, they release germs into the air that may infect others nearby. Face masks are part of an infection control strategy to eliminate cross-contamination.

How to Wear Mask



- Before putting on the mask, wash your hands for at least 20 seconds with soap and water, or rub your hands together thoroughly with hand sanitizer.
- The nose clip part should be facing up and positioned against bridge of your nose
- Hold the mask by both ear loops and place one loop over each ear.
- Mold the bendable nose clip to the shape of your nose by pinching and pressing down on it with your fingers.
- Pull the bottom of the mask over your mouth and chin.
- Be sure the mask fits snugly.
- Do not touch the mask once in position.
- If the mask gets soiled or damp, replace it with a new one.

**Images for representation purpose only. Actual product may be slightly different.*